Stop worrying about your cells overheating



AutoMate Scientific's new ThermoClamp®-3 temperature controller can maintain the temperature of your cells in electrophysiology and imaging experiments by heating your solutions, the perfusion chamber, or both.



- Combination inline heater plus multi-channel focal drug delivery
- Advanced auto-tuning temperature lock with fuzzy logic
- Ultra low-noise for physiology research
- Many inputs and outputs for computer control and recording
- Single and dual controllers available
- New low price starting at \$999



Thermal Runaway means a temperature controller has lost contact with the feedback sensor and keeps increasing the heat until a maximum temperature is reached because the controller thinks your cells are actually too cold.


The ThermoClamp can be used with one or two temperature sensors. There is always a "Safety" sensor on the heater, and an optional "Control" sensor in your bath. Using two sensors is more accurate, but causes thermal runaway if the perfusion stops or the bath sensor comes out of the chamber.


The patent-pending ThermoClamp-3 is the world's first and only physiology temperature controller that detects and prevents thermal runaway.

Other temperature controllers

AutoMate Scientific ThermoClamp®-3 controller

Thermal Runaway when the bath sensor falls out or you change coverslips or dishes:






Thermal Runaway when perfusion flow stops:




PRODUCTS GIMBH

for Research in Life Sciences



SCIENCE



Science Products GmbH

Hofheimer Straße 63 D- 65719 Hofheim Phone: +49 6192 90 13 96

Fax: +49 6192 90 13 98

+43 7207 7 54 55

Email: info@science-products.com

+48 12 363 80 52

in collaboration with:

Science Products Trading SPT AG

Käferholzstr. 142 CH- 4058 Basel

Phone +41 43 488 05 61 Fax: +41 43 488 05 62

